


Nálevka dělicí kónická, s PTFE ventilem Rotaflo[®], bez zátky

Objem ml	Zábrus	Průměr výpusti mm	Kat. číslo	Kč za kus
100	14/23	3,0	1415.9280	1 414,-
100	19/26	3,0	1415.9282	1 387,-
250	24/29	3,0	1415.9284	1 794,-
500	24/29	6,0	1415.9286	1 862,-



01

02

03

Nálevka dělicí kuželová podle Squibba, s NZ zátkou a teflonovým NZ kohoutem

Objem ml	SJ	d mm	l mm	Počet ks v balení	Kat. číslo	Kč za kus
50	19/26	9	70	4*	1415.7030	655,-
100	19/26	9	70	4*	1415.7032	692,-
250	29/32	9	70	2*	1415.7034	729,-
500	29/32	10	70	2*	1415.7036	867,-
1000	29/32	10	70	2*	1415.7038	1 507,-
2000	29/32	13	70	4*	1415.7040	1 280,-

* lze objednávat i po jednotlivých kusech

SIMAX


04

05

06

Nálevka dělicí podle Squibba, hnědé sklo, se skleněným kohoutem a plastovou zátkou

Objem ml	Zábrus	Průměr výpusti mm	Počet kusů v balení	Kat. číslo	Kč za balení
50	19/26	2,5	2	1415.9040	2 688,-
100	19/26	2,5	2	1415.9042	3 054,-
250	29/32	4,0	2	1415.9044	3 760,-
500	29/32	4,0	2	1415.9046	4 356,-
1000	29/32	6,0	2	1415.9048	5 740,-
2000	29/32	6,0	1	1415.9050	3 380,-



07

08

09

Nálevka dělicí podle Squibba, hnědé sklo, s PTFE kohoutem a plastovou zátkou

Objem ml	Zábrus	Průměr výpusti mm	Počet kusů v balení	Kat. číslo	Kč za balení
50	19/26	2,5	2	1415.9060	3 200,-
100	19/26	2,5	2	1415.9062	3 252,-
250	29/32	4,0	2	1415.9064	3 896,-
500	29/32	4,0	2	1415.9066	4 684,-
1000	29/32	6,0	2	1415.9068	5 794,-
2000	29/32	6,0	1	1415.9070	3 382,-



10

11

12