


Baňka odměrná s vyhnutým okrajem, bez zátky

Fisherbrand[®] SIMAK

třída A, šaržový certifikát

Objem ml	Max. odchylka ± ml	d mm	h mm	Počet ks v balení	Kat. číslo
5	0,025	22	70	10*	1351.1018
10	0,025	27	90	10*	1351.1019
20	0,04	40	110	5*	1351.5022
25	0,04	40	110	10*	1351.1023
50	0,06	50	140	10*	1351.1025
100	0,10	60	170	10*	1351.1030
200	0,15	75	210	10*	1351.1037
250	0,15	80	220	10*	1351.1038
500	0,25	100	260	10*	1351.1043
1000	0,40	125	300	10*	1351.1044
2000	0,60	160	370	6*	1351.1046

* lze objednávat i po jednotlivých kusech



Baňka odměrná s vyhnutým okrajem, bez zátky


SIMAK

třída K, jde o metrologicky ověřené baňky třídy A s ověřovacím listem.

Baňka je označena metrologickou značkou.

Objem ml	Max. odchylka ± ml	d mm	h mm	Počet ks v balení	Kat. číslo
5	0,025	22	70	10*	1351.2018
10	0,025	27	90	10*	1351.2019
25	0,04	40	110	10*	1351.2023
50	0,06	50	140	10*	1351.2025
100	0,10	60	170	10*	1351.2030
200	0,15	75	210	10*	1351.2037
250	0,15	80	220	10*	1351.2038
500	0,25	100	260	10*	1351.2043
1000	0,40	125	300	10*	1351.2044
2000	0,60	160	370	6*	1351.2046

* lze objednávat i po jednotlivých kusech




Baňka odměrná se zábrusem a skleněnou dutou zátkou

Fisherbrand[®] SIMAK

třída A, šaržový certifikát

Objem ml	Max. odchylka ± ml	d mm	h mm	SJ mm	Počet ks v balení	Kat. číslo
5	0,025	22	70	7/16	5*	1352.1018
10	0,025	27	90	10/19	5*	1352.1019
25	0,04	40	110	10/19	5*	1352.1023
50	0,06	50	140	12/21	5*	1352.1025
100	0,10	60	170	12/21	5*	1352.1030
200	0,15	75	210	14/23	2*	1352.1037
250	0,15	80	220	14/23	2*	1352.1038
500	0,25	100	260	19/26	2*	1352.1043
1000	0,40	125	300	24/29	2*	1352.1044
2000	0,60	160	370	29/32	6*	1352.1046

* lze objednávat i po jednotlivých kusech



Baňka odměrná se zábrusem a skleněnou zátkou

SIMAX

třída K, jde o metrologicky ověřené baňky třídy A s ověřovacím listem

Baňka je označena metrologickou značkou.



Objem ml	Max. odchylka ± ml	d mm	h mm	SJ mm	Počet ks v balení	Kat. číslo
5	0.025	22	70	7/16	10*	1352.3018
10	0.025	27	90	7/16	10*	1352.3019
25	0.04	40	110	10/19	10*	1352.3023
50	0.06	50	140	12/21	10*	1352.3025
100	0.10	60	170	12/21	10*	1352.3030
200	0.15	75	210	14/23	10*	1352.3037
250	0.15	80	220	14/23	10*	1352.3038
500	0.25	100	260	19/26	10*	1352.3043
1000	0.40	125	300	24/29	10*	1352.3044
2000	0.60	160	370	29/32	6*	1352.3046


* lze objednávat i po jednotlivých kusech

Baňka odměrná se zábrusem a plastovou zátkou



Fisherbrand SIMAX

třída A, šaržový certifikát



Objem ml	Max. odchylka ± ml	d mm	h mm	SJ mm	Počet ks v balení	Kat. číslo
5	0.025	22	70	10/19	5*	1353.1018
10	0.025	27	90	10/19	5*	1353.1019
20	0.04	40	110	10/19	5*	1353.1022
25	0.04	40	110	10/19	5*	1353.1023
50	0.06	50	140	12/21	5*	1353.1025
100	0.10	60	170	12/21	5*	1353.1030
200	0.15	75	210	14/23	2*	1353.1037
250	0.15	80	220	14/23	2*	1353.1038
500	0.25	100	260	19/26	2*	1353.1043
1000	0.40	125	300	24/29	2*	1353.1044
2000	0.60	160	370	29/32	6*	1353.1046
5000	1.80	-	-	34/35	3*	1353.1049
10000	2.00	-	-	45/90	1*	1353.1053


* lze objednávat i po jednotlivých kusech

Baňka odměrná se zábrusem a plastovou zátkou

SIMAX

třída K, jde o metrologicky ověřené baňky třídy A s ověřovacím listem.

Baňka je označena metrologickou značkou.



Objem ml	Max. odchylka ± ml	d mm	h mm	SJ mm	Počet ks v balení	Kat. číslo
5	0.025	22	70	7/16	10*	1353.3018
10	0.025	27	90	7/16	10*	1353.3019
25	0.04	40	110	10/19	10*	1353.3023
50	0.06	50	140	12/21	10*	1353.3025
100	0.10	60	170	12/21	10*	1353.3030
200	0.15	75	210	14/23	10*	1353.3037
250	0.15	80	220	14/23	10*	1353.3038
500	0.25	100	260	19/26	10*	1353.3043
1000	0.40	125	300	24/29	10*	1353.3044
2000	0.60	160	370	29/32	6*	1353.3046

* lze objednávat i po jednotlivých kusech

Baňka odměrná se zábrusem a skleněnou zátkou

BLAUBRAND® třída A, šaržový certifikát

Objem ml	Max. odchylka ± ml	d mm	h mm	SJ mm	Počet kusů v balení	Kat. číslo
5W	0,04	22	70	10/19	2	1353.4250
10W	0,04	27	90	10/19	2	1353.4252
20	0,04	39	110	10/19	2	1353.4254
20W	0,06	39	105	12/21	2	1353.4256
25	0,04	40	110	10/19	2	1353.4258
25W	0,06	40	110	12/21	2	1353.4260
50	0,06	50	140	12/21	2	1353.4262
50W	0,1	50	140	14/23	2	1353.4264
100	0,1	60	170	12/21	2	1353.4266
100	0,1	60	170	14/23	2	1353.4268
200	0,15	75	210	14/23	2	1353.4270
250	0,15	80	220	14/23	2	1353.4272
500	0,25	100	260	19/26	2	1353.4274
1000	0,4	125	300	24/29	1	1353.4276
2000	0,6	160	370	29/32	1	1353.4278
5000	1,2	215	475	34/35	1	1353.4280

W - širokohlávké provedení



Baňka odměrná se zábrusem a plastovou zátkou

BLAUBRAND® třída A, šaržový certifikát

Objem ml	Max. odchylka ± ml	d mm	h mm	SJ mm	Počet kusů v balení	Kat. číslo
5W	0,04	22	70	10/19	2	1353.4200
10W	0,04	27	90	10/19	2	1353.4202
20	0,04	39	110	10/19	2	1353.4204
20W	0,06	39	105	12/21	2	1353.4206
25	0,04	40	110	10/19	2	1353.4208
25W	0,06	40	110	12/21	2	1353.4210
50	0,06	50	140	12/21	2	1353.4212
50W	0,1	50	140	14/23	2	1353.4214
100	0,1	60	170	12/21	2	1353.4216
100	0,1	60	170	14/23	2	1353.4218
200	0,15	75	210	14/23	2	1353.4220
250	0,15	80	220	14/23	2	1353.4222
500	0,25	100	260	19/26	2	1353.4224
1000	0,4	125	300	24/29	1	1353.4226
1000W	0,6	125	300	29/32	1	1353.4228
2000	0,6	160	370	29/32	1	1353.4230
5000	1,2	215	475	34/35	1	1353.4232
10000	2	270	570	45/40	1	1353.4234

W - širokohlávké provedení



Baňka odměrná se zábrusem a plastovou zátkou

BLAUBRAND® ETERNA třída A, šaržový certifikát

Objem ml	Max. odchylka ± ml	d mm	h mm	SJ mm	Počet kusů v balení	Kat. číslo
5W	0,04	22	70	10/19	2	1353.4300
10W	0,04	27	90	10/19	2	1353.4302
20	0,04	39	110	10/19	2	1353.4304
20W	0,06	39	105	12/21	2	1353.4306
25	0,04	40	110	10/19	2	1353.4308
25W	0,06	40	110	12/21	2	1353.4310
50	0,06	50	140	12/21	2	1353.4312
50W	0,1	50	140	14/23	2	1353.4314
100	0,1	60	170	14/23	2	1353.4316
200	0,15	75	210	14/23	2	1353.4318
250	0,15	80	220	14/23	2	1353.4320
500	0,25	100	260	19/26	2	1353.4322
1000	0,4	125	300	24/29	1	1353.4324
2000	0,6	160	370	29/32	1	1353.4326

W - širokohlávké provedení



Baňka odměrná se zábrusem a plastovou zátkou

HIRSCHMANN®



třída A, hnědá graduace, včetně šaržového certifikátu dle DIN EN ISO 1042

Objem ml	Max. odchylka ± ml	h ± 5 mm	SJ mm	Počet ks v balení	Kat. číslo
5	0,025	70	7/16	2	1353.6018
10	0,025	90	7/16	2	1353.6019
20	0,04	110	10/19	2	1353.6022
25	0,04	110	10/19	2	1353.6023
50	0,06	140	12/21	2	1353.6025
100	0,10	170	12/21	2	1353.6030
200	0,15	210	14/23	2	1353.6037
250	0,15	220	14/23	2	1353.6038
500	0,25	260	19/26	2	1353.6043
1000	0,40	300	24/29	1	1353.6044
2000	0,60	370	29/32	1	1353.6046
5000	1,20	470	34/35	1	1353.6049
10000	2,00	570	45/40	1	1353.6053

Baňka odměrná se zábrusem a skleněnou dutou zátkou

HIRSCHMANN®




třída A, modrá graduace, včetně šaržového certifikátu dle DIN EN ISO 1042

Objem ml	Max. odchylka ± ml	h ± 5 mm	SJ mm	Počet ks v balení	Kat. číslo
5	0,025	70	7/16	2	1353.6118
10	0,025	90	7/16	2	1353.6119
20	0,04	110	10/19	2	1353.6122
25	0,04	110	10/19	2	1353.6123
50	0,06	140	12/21	2	1353.6125
100	0,10	170	12/21	2	1353.6130
200	0,15	210	14/23	2	1353.6137
250	0,15	220	14/23	2	1353.6138
500	0,25	260	19/26	2	1353.6143
1000	0,40	300	24/29	1	1353.6144
2000	0,60	370	29/32	1	1353.6146
5000	1,20	470	34/35	1	1353.6149
10000	2,00	570	45/40	1	1353.6153

Baňka odměrná se zábrusem a plastovou zátkou, hnědé sklo

HIRSCHMANN®



třída A, hnědé borosilikátové sklo

Objem ml	Max. odchylka ± ml	d mm	h mm	SJ mm	Počet ks v balení	Kat. číslo
5	0,025	22	70	7/16	2	1353.5011
10	0,025	27	90	7/16	2	1353.5012
25	0,04	40	110	10/19	2	1353.5013
50	0,06	50	140	12/21	2	1353.5014
100	0,1	60	170	12/21	2	1353.5015
200	0,15	75	210	14/23	2	1353.5016
250	0,15	80	220	14/23	2	1353.5017
500	0,25	100	260	19/26	2	1353.5018
1000	0,4	125	300	24/29	1	1353.5019
2000	0,6	160	370	29/32	1	1353.5020

Baňka odměrná potažená plastem, se zábrusem a plastovou zátkou

BLAUBRAND® PURprotect třída A, šaržový certifikát

Objem ml	Max. odchylka ± ml	d mm	h mm	SJ mm	Počet kusů v balení	Kat. číslo
50W	0,1	50	140	14/23	2	1353.4400
100	0,1	60	170	14/23	2	1353.4402
200	0,15	75	210	14/23	2	1353.4404
250	0,15	80	220	14/23	2	1353.4406
500	0,25	100	260	19/26	2	1353.4408
1000	0,4	125	300	24/29	1	1353.4410

W - širokohrdlé provedení



Baňka odměrná trapézoidní, se zábrusem a plastovou zátkou

BLAUBRAND® třída A, šaržový certifikát

Objem ml	Max. odchylka ± ml	d mm	h mm	SJ mm	Počet kusů v balení	Kat. číslo
1	0,025	13	65	7/16	2	1353.4000
2	0,025	17	70	7/16	2	1353.4002
5	0,025	22	70	7/16	2	1353.4004
10W	0,04	27	90	10/19	2	1353.4006
20	0,04	39	110	10/19	2	1353.4008
25	0,04	40	110	10/19	2	1353.4010
50	0,06	50	140	12/21	2	1353.4012

W - širokohrdlé provedení

