

XS pH 60 VioLab benchtop pH meter - Electrode 201 T

Code: 50101092



TECHNICAL DATA SHEET

Supplied as standard

Supplied as standard

Electrode holder stand, coloured buffer solutions, USB cable, multi-socket power supply and instruction manual.

Performance

Accuracy (pH)	± 0,02 pH pH
Temperature accuracy	± 0,5 °C °C
Temperature compensation	0,0 ... +100,0 °C °C
Measuring range mV	-1000 ... +1900 mV mV
pH measuring range	-2,00 ... +16,00 pH pH
temperature measurement range	-10,0 ... +110,0 °C °C
pH resolution	0,1 / 0,01 pH pH
Temperature Resolution	0,1 °C °C
Working temperature	0 ... +45 °C °C

Features

PC Communication	Micro USB
Printer connection	RS 232
Criteria for stability of the measure	Low - Med - High
Display	High-definition, backlit colour LCD
AQA	Yes
IP rating	54 IP IP
Display	Yes with icon
Inputs	BNC - RCA
Memory / Datalogger	1000
PH calibration points	1 ... 3 automatic values - 2 user values
Calibration Report	Yes, printable
Recognised pH buffers	USA - NIST

Via della Meccanica, 25, 41012 Carpi (MO) - Italy

+39 059 653274 - info@giorgiobormac.com

Display

Si

Feeds

Power supply

100 - 240 V V

Dimensions and weight

Dimensions (WxDxH) mm

160 x 140 x 45 mm

Weight

0.38 kg Kg
